



Alcohol-based rapid disinfectant for the disinfection of non-invasive medical products.

terralin® liquid

Our Plus

- aldehyde free
- ready-to-use
- rapid and effective in many applications
- dries rapidly, leaving no streaks

Application areas

Alcohol-based rapid disinfection of medical devices in all areas and wipeable surfaces with an increased risk of infection and where short contact times are required, e.g.:

- Patient nearby surfaces
- Patient treatment units
- Examination couches
- Operating tables and nearby work surfaces
- Surfaces of medical equipment

Instructions for use

Apply to surfaces undiluted, wipe and allow to take effect. Make sure to wet surfaces completely and keep them wet for the whole exposure time. Maximum application rate 50 ml/sq.m. Make sure that all visible soiling is removed before disinfection and the lid is tightly closed after use. Do not treat surfaces that are sensitive to alcohol, e.g. acrylic glass. Not suitable for the disinfection of invasive medical devices.

Microbiological efficacy

| Efficacy | Concentration | Contact time |
|--|---------------|--------------|
| bactericidal <i>in accordance with VAH Guidelines with short contact times</i> | ready-to-use | 1 min. |
| MRSA | ready-to-use | 1 min. |
| tuberculocidal <i>EN14348</i> | ready-to-use | 1 min. |

| Efficacy | Concentration | Contact time |
|--|---------------|--------------|
| enveloped viruses (incl. HIV, HBV, HCV) <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i> | ready-to-use | 30 sec. |
| norovirus <i>EN14476</i> | ready-to-use | 1 min. |
| polyoma SV40 <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i> | ready-to-use | 10 min. |
| rotavirus | ready-to-use | 30 sec. |
| adenovirus (type 5) <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i> | ready-to-use | 5 min. |

Certificates

- IHO listed
- VAH certificate



Product data

Composition:

100 g contain the following active ingredients:
25 g Ethanol (94 %), 35 g Propan-1-ol

Chemical-physical data

| | |
|--------------------|------------------------------------|
| Color | colourless |
| Density | ca. 0,89 g/cm ³ / 20 °C |
| Flash point | 27 °C / Method : DIN 51755 Part 1 |
| Form | liquid |
| pH | Not applicable |
| Viscosity, dynamic | not determined |

Special advice

Always read the label and product information before use.

The product is highly compatible with metals and synthetics (except acrylic glass and alcohol-sensitive varnishes). Please follow carefully any national fire and explosion protection guidelines when using alcohol-based disinfectants.

Information for order

| Item | Delivery form | Item no. |
|---|---------------|------------|
| terralin [®] liquid 1 l bottle | 10/Carton | on request |
| terralin [®] liquid 10 l | 1/Canister | on request |

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Accessories

| Accessories | Item no. |
|---|----------|
| 5 l/10 l-Canister key | 135810 |
| 5 l/10 l-Canister tap | 135501 |
| Canister dosing feeder 5/10l (20ml / pump stroke) | 117101 |
| Compressed-air nozzle | 120001 |
| Decanter funnel | 117901 |
| trigger-sprayer for 500 ml /1000 ml | 180124 |

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Care

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1 AT
Phone +44 (0) 1142 - 5435 - 00
Fax +44 (0) 1142 - 5435 - 01
www.schuelke.co.uk
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H.
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 5232501-0
Fax +43 (0) 1 - 5232501-60
www.schuelke.at
office.austria@schuelke.com